

# Study List (Studies Screen)

Studies received from Tricefy are populated on your Studies screen. Select **Studies** at the top of the page to access this screen.

## Screen Layout

The information on your Studies Screen is determined by the column headers. By default, these include:

- **Thumbnails:** ► indicates the thumbnails are collapsed. Select the image icon in the column header to toggle thumbnails or can open thumbnails for individual studies by selecting the ► buttons
- **Patient's full name and study date:** This data is received from the ultrasound system
- **Status:** These are created by your Account Administrator
  - To change the status of a study, double-click the status
- **Study Details:** Indicates the type of files included in the study, such as:
  - 🖼️ indicates the number of images
  - 🎞️ indicates the cine clips
  - 👤 indicates the number of collaborators
  - 📄 indicates the number of consults (shares with doctors)
  - 📝 indicates the number of notes, including annotations
  - 📧 indicates how many times a link was sent to the patient

Column header can reorganize your studies. Select a column header to group studies together and list them in alphanumeric or descending/ascending order.

For example:

- Selecting the Full Name column header will reorganize the studies, so they are listed in alphabetical order by name (A - Z)
- Selecting Full Name again will reverse the order (Z - A)
- Selecting the Status column header will list the studies grouped together by status

## Page Settings

Columns are customizable. You can change, remove, and add additional column headers, such as the sonographer, patient's date of birth, procedure, etc. by selecting the **settings cog** at the top of the Studies screen.

- Check **Display Thumbnails** to always show thumbnails for each study
- Check **Thumbnails in Reverse Order** to always display thumbnails in reverse order (with reports first, if applicable)
- Uncheck **Study List Scrolls Horizontally** to disable scrolling (right to left) through thumbnails for individual studies
- Select or deselect options to enable or disable column headers
  - A **green** check indicates the column is visible

## Reordering Column Headers

To change the order of your column headers, drag-and-drop header listings in a different order within the Page Settings.

## Column Width

By default, the column widths fit to the necessary size. To change the width:

1. Hover over the line that separates two columns until it turns into a double-sided arrow
2. Hold while slightly moving the arrow to make the columns wider or smaller

# Basic Operation

When you select a study (examination) on your Study List, that study will turn **blue**.

If the thumbnails are expanded, they are all selected (you will know they are selected because they will be outlined in **green**). This means that when you send them to a patient or share them with a doctor, all images are sent.

When a study is selected, additional options appear at the top of your screen:

- **Open in viewer:** Opens the examination in the Viewer, where you can zoom, view multiple images simultaneously, play video clips, and make notes
- **Send to:** Provides options for sending the examination to a patient, doctor (such as a referring physician) or external server, like an EMR
- **Action menu:** Displays options for attaching, moving, or deleting the examination, as explained below
- **Label icon:** Displays labels that can be applied to the examination

If Tricefy's adaptive reporting is enabled, there is also an option to create a report.

## Action Menu

The Action Menu (☰) is at the top of the studies screen, consults screen and Viewer. It is only visible when a study is selected:

- **Attach Files** \*: Select a file to attach to the study
- **Edit Study** \*: Assign the study to a different patient (if their name was misspelled in the ultrasound system)
- **Move:** Transfer or copy a study to a different Tricefy account (only available for members of multiple accounts)
- **Download Studies:** Download a study from Tricefy to your computer
- **Delete:** Remove a study from Tricefy. You can resend the study using your imaging system

- **Clear Selection:** De-selects the study/studies
- **Select All:** Selects all the thumbnails in a study; only available when using the Viewer

\* Not available when multiple studies are selected

## Send to Menu

The **Send To** menu is located at the top of your Viewer and Studies screen (if a study is selected)

- **Send to Patient:** Sends an examination to a patient
  - They will receive a link to their own special webpage hosting their images
- **Send to Doctor:** Sends an examination to a doctor, such as a referring physician
  - They will receive a link to view and comment on the exam.
- **Send to DICOM Server:** Sends an examination to an external server, such as a PACS
  - Tricefy integrates with outside software and servers to accommodate efficient workflows

On the Consults screen, the **Send to DICOM Server** is the only option in this menu.

- To send the exam to patient, return to your Studies screen by selecting **Studies** from the navigational options
- To send the exam to additional doctors, use the side information panel on the Consults screen

## Side Information Panel

You can attach files, apply labels, and view study information using the Side Information Panel.

If the side panel is opened on the Studies screen and a study is not selected, the panel displays a time-stamped history log of your clinic's activity.

- **Account Activity tab:** Displays a history log that includes all members of your clinic
- **My Activity tab:** Displays a history log of only your actions

Once a study is selected, the information in the panel will display information specific to that study.

This information is organized into six tabs:

### Info Tab

The **Info tab** shows the date of the study, when the study expires, patient

information (name, gender, date of birth), procedure information (physicians, procedure) and the size of the study.

### Attaching Files (side panel – Info Tab)

To attach files, such as a report, PDF, or image file, select **Choose Files** and select files from your computer/network.

### Reports Tab

The **Reports tab** (available only for our reporting customers) lists pending and finalized reports. Select a pending report to continue it; select a finalized report to view or update it.

### Labels Tab

Labels help filter, sort, and search for studies. To add (or remove) a label for the study, select the **+** button. Previously used labels are pre-populated; check or uncheck as applicable. To add a new label, type into the text field, followed by **Enter** on the keyboard.

### Shares Tab

The **Shares tab** displays a log of when and where the study was previously shared.

- If the study was shared with a patient, it will be logged under **Patient Links**
  - If the patient accidentally deletes the images or needs them sent again resend them by selecting the drop-down arrow. This method will send the same link to the same contact information.
  - If the link expired or the images need to be sent using different contact information, share the images with the patient again using the **Send to Patient** feature.

If the study was sent to a colleague/another doctor, it will be logged under **Consults**. It will let you know when your colleague viewed the study. If they have yet to view the study, it will say, **Unread**.

### Notes Tab

The **Notes tab** displays time-stamped notes and a text box for creating new notes. Type a comment into the text box, followed by **Add Note**. This note will be visible to the other members of your clinic. Select the **X** to erase your note.

### Activity Tab

The **Activity tab** displays a time-stamped log for actions performed by your clinic.