

GE Healthcare

VividE80-VividE90-VividE95  
VividS60-VividS70  
EchoPAC  
Quick Reference Guide  
Connectivity



# Vivid E80/E90/E95/S60/S70 Quick Ref. Guide Index

## Connectivity - Basics

- Connections
- Getting started
- Dataflows
- Compatibility
- Definitions
- Network Settings
- DICOM Activities

## Connectivity - Configure & Connect to .....

- Worklist
- DICOM Storage
- DICOM CD/DVD
- Query Retrieve
- XML Export
- Insite ExC
- DICOM Printer
- Network Printer
- DICOM USB
- Remote Archive
- DICOM Spooler
- Backup & Restore
- Disk Management
- EchoPAC PC/SWO
- MPEGVue Export
- Viewpoint

## Connectivity – Quick Aid

- Checklist
- Troubleshooting
- Repair DICOMDIR
- Local Archive Corrupt

# Connections

## Vivid E80/E90/E95



Vivid E80/E90/E95 rear view




**Ethernet**  
Network  
1000 Base-TX  
Ethernet IEEE 802.3

USB



Printer

USB 2.0 only (not 3.0)



External monitor  
Connector Type :  
**DVI-D**

**Note!** The External Monitor output **signal** is DVI-D (Digital),  
Do **not** connect an analog monitor to this port

# Connections

## Vivid S60/S70

### Insulated USB



USB 1.0  
Devices w/  
external power  
Printer

### Wifi (Optional)



USB 2.0 only  
(not 3.0) . Lower  
Socket

WiFi Network Adaptor  
(Standards: 802.11b, 802.11g,  
802.11n)

### USB Dual Socket




USB 2.0 only  
(not 3.0) . Devices  
w/o external power  
Memory devices

### Ethernet



Network  
1000 Base-TX  
Ethernet IEEE 802.3



External monitor  
Connector Type :  
**DVI-D**



Vivid S60/S70 rear view


**Note!** The External Monitor output **signal** is DVI-D (Digital),  
Do **not** connect an analog monitor to this port

# Connections

## Network


### Direct Cable Connection

- Connect directly with a Crossover cable




### Peer-to-Peer Network

- Connect to HUB/Switch with a network cable



### Hospital Network

- Connect to wall outlet with a network cable



### Network requirements

- To achieve a successful connection to the different devices, the following is required:

#### DICOM device

- IP Address
- AE title
- Portnumber
- Image Settings
- Dicom SR Settings

#### Remote Device

- IP Address
- Host Name

#### Network Printer

- IP Address
- Host Name

- Devices in the same network must have the 3 first segments identical, the 4th must be unique

Typical IP Addr :   
Identical Unique  
 aaa.bbb.ccc.xxx

# Definitions

**DICOM** - Digital Imaging and Communications in Medicine

## Terminology

**Application Entity Title** – the externally known name of an *Application Entity*, used to identify a DICOM application to other DICOM applications on the network.

**Application Entity (AE)** – an end point of a DICOM information exchange, including the DICOM network or media interface software; i.e., the software that sends or receives DICOM information objects or messages. A single device may have multiple Application Entities.

**Joint Photographic Experts Group (JPEG)** – a set of standardized image compression techniques, available for use by DICOM applications.

**Service Class Provider (SCP)** – role of an *Application Entity* that provides a DICOM network service; typically, a server that performs operations requested by another *Application Entity* (*Service Class User*). Examples: Picture Archiving and Communication System (image storage SCP, and image query/retrieve SCP), Radiology Information System (modality worklist SCP).

**Service Class User (SCU)** – role of an *Application Entity* that uses a DICOM network service; typically, a client. Examples: imaging modality (image storage SCU, and modality worklist SCU), imaging workstation (image query/retrieve SCU)

**Service/Object Pair (SOP) Class** – the specification of the network or media transfer (service) of a particular type of data (object); the fundamental unit of DICOM interoperability specification. Examples: Ultrasound Image Storage Service, Basic Grayscale Print Management.

**Tag** – a 32-bit identifier for a data element, represented as a pair of four digit hexadecimal numbers, the “group” and the “element”. If the “group” number is odd, the tag is for a private (manufacturer-specific) data element. Examples: (0010,0020) [Patient ID], (07FE,0010) [Pixel Data], (0019,0210) [private data element]

**Unique Identifier (UID)** – a globally unique “dotted decimal” string that identifies a specific object or a class of objects; an ISO-8824 Object Identifier. Examples: Study Instance UID, SOP Class UID, SOP Instance UID.

# Definitions

## Terminology

**DICOM** - Digital Imaging and Communications  
in Medicine

DNS	Domain Name System
FSC	File-Set Creator
FSU	File-Set Updater
FSR	File-Set Reader
HIS	Hospital Information System
HL7	Health Level 7 Standard
IHE	Integrating the Healthcare Enterprise
IOD	Information Object Definition
IPv4	Internet Protocol version 4
IPv6	Internet Protocol version 6
ISO	International Organization for Standards
JPEG	Joint Photographic Experts Group
LUT	Look-up Table
MPPS	Modality Performed Procedure Step
MSPS	Modality Scheduled Procedure Step
MTU	Maximum Transmission Unit (IP)
MWL	Modality Worklist
RDCM	Raw Data Conversion Module

# Definitions

**DICOM** -Digital Imaging and Communications  
in Medicine

## Terminology

### **DICOM Storage Commit - what is it ?**

DICOM Storage Commitment is a way for an Image Storage Service to inform your scanner that it has successfully stored your images without corruption. With this confirmation, you may decide to delete the images from your scanner.

### **DICOM MPPS - what is it ?**

A DICOM Performed Procedure, also known as Modality Performed Procedure Step (MPPS) is a way to inform your billing system that an ultrasound procedure has been completed and can be billed.

A typical setup has your Worklist Service handle MPPS communications so that, for example, you can retrieve today's scheduled exams from, and report completed exams to, the same system. In basics, it sends confirmation to the MWL when the exam has started/finished

### **DICOM SR - what is it ?**

Structured Reporting provides for transmission of exam measurements and calculations in a separate file. It is always associated with a DICOM Image Storage service. The receiver must be able to interpret and display the data in the Structured Report.

### **DICOM SR Private Tags- what is it ?**


Within DICOM, manufacturers have created their own tags, containing specific data elements for their own use. This is called a private tag, and will generally not be recognized by another application reading their DICOM file. Private tags are useful when manufacturers want to convey information not defined within the DICOM standard, to be read generally by another instance of their application, or another application within their design portfolio




# Getting Started

## Go to Connectivity Page

1. Switch on the circuit breaker on the rear of the unit, press the on/off button to start



Vivid E80/E90/E95 Circuit Breaker location




Vivid S60/S70 Circuit Breaker location

Wait for the Application to start and complete

2. Select  located in upper right corner on the Touchscreen to enter **Utility** menu


3. Select  or press «F2» on the Alphanumeric keyboard

A login window will appear :



Log in window

If purpose is to change Dataflows and TCP/IP, log in as **ADM**, Administrator, with appropriate password



Touch Screen

# Getting Started


## Go to Connectivity Page

4. Select  button to get to Connectivity Page

Printer, DICOM Printer,  
Output buttons  
Recorder set up

Dataflow/workflow  
selections

Dataflow/workflow  
Overview,  
Input & Output



TCP/IP settings

Store Image directly or  
to buffer

Hide /show dataflow

Connectivity Page

# Network Settings

## TCP/IP

Select **Tcpip** tab to get to TCP/IP Page

### My Computer

Computer name, IP Address from Windows settings.

AE Title & Port No according to hospital IT policy

### Remote Path

Use UNC syntax for shared folder:

[\\SERVER\SharedFolder](#)

### Configurable Remote Path User

If the default user cannot be used, insert User and Password to connect to the Remote Path. User can be 'DomainName\User' syntax if the shared folder is on a Network Domain

### Save Settings

Always select this button whenever there is a change made on this page.

### TCP/IP settings

The screenshot shows the 'CONNECTIVITY' window with the 'Tcpip' tab selected. Red boxes and arrows highlight specific areas: 'My Computer' fields (Computer Name, IP-Address, AE Title, Port No), 'Server Config' (Servers list), 'Remote Path' (Remote Path field, Check button), 'Configurable Remote Path User' (User/Password fields, Note), 'Save settings' button, 'Network Settings' button, and 'DICOM' section (Detailed DICOM Log checkbox). Red lines also connect the explanatory text blocks to their corresponding fields in the interface.

Section	Field/Control	Value	
My Computer	Computer Name	HCE-HNYKYW1.logon.ds.ge.com	
	IP-Address	192.168.1.16	
	AE Title	HCE-HNYKYW1	
	Port No	104	
Server Config	Servers	ECHOPAC-000000 ) 10.0.0.4 DICOMSERVER ) 10.0.0.5 HL7 ) 10.0.0.7 NAS ) 192.168.1.60	
	Buttons	Modify, Add, Remove	
	Remote Path	Remote Path	\\SERVER\SHARE
	Buttons	Check	
Configurable Remote Path User	User		
	Password		
DICOM	Detailed DICOM Log	<input type="checkbox"/>	

### Server Config

Add, Modify and Remove Servers that will be used in the different Dataflows

### Check

Network share connectivity test with default user and password or «Configurable Remote Path User» as secondary log in credential

### Network Settings

Select this button to open Windows Network setup to change the IP address on the Scanner

### Detailed DICOM Log

Enable when troubleshooting on DICOM devices only

TCP IP Page

# Dataflow

## Concept

Communication between the Vivid E80 / E90 / E95 ultrasound Scanner/EchoPAC and other information providers on the network takes the form of dataflows. Each dataflow defines the transfer of patient information and images from an input source to the unit, and from the unit to one or several output sources.

A dataflow is a set of pre-configured settings. Selecting a dataflow will automatically customize the unit to work according to the settings associated with this dataflow.

Select the  tab on the Connectivity Page to get to the Dataflow page

# Dataflows

## Selections

### LocalArchive-Int.HD

Internal HD (Vivid)




The local database is used for patient archiving. Images are stored to the internal hard drive.

### LocalArchive - Int HD/DICOM Server

Internal HD (Vivid)

DICOM Server




The local archive is used for patient archiving. Images are stored to the internal hard drive and to a DICOM server.

### Worklist/LocalArchive-DICOMServer/Int.HD

DICOM Worklist

Internal HD (Vivid)

DICOM Server



Search in the DICOM Modality Worklist, the patient found is copied into local database. The patient information and the examination results are stored to the local the database. Images are stored to a DICOM Server and to an image volume on the internal hard drive. Some of the measurements are stored if DICOM SR is turned on


### RemoteArch-RemoteHD

Vivid

REMOTE ARCHIVE

EPSWO  
ImageVault

REMOTE ARCHIVES




A remote database is used for patient archiving. Images are stored to a network image volume

### Remote Archive - Remote HD/DICOM Server

REMOTE ARCHIVE

DICOM Server



A remote database is used for patient archiving. Images are stored to a network image volume and to a DICOM server. Some of the measurements are stored if DICOM SR is turned on


### Worklist/RemoteArchive-DICOMServer/RemoteHD

DICOM Worklist

Vivid

REMOTE ARCHIVE

DICOM Server




Search in a DICOM Modality Worklist, the patient found is copied into a remote database. The patient information and the examination results are stored to a remote database. Images are stored to a DICOM server and to an image network volume. Some of the measurements are stored if DICOM SR is turned on.

# Dataflows

## Selections

### Worklist/LocalArchive-Int.HD




Search in the DICOM Modality Worklist, the patient found is copied into the local database. The patient information and examination results are stored to the local database. Images are stored to the internal hard drive..

### LocalArchive-Int.HD/eVue




The local database is used for patient archiving. Images are stored to internal hard drive and an MPEG exam is created in the configured destination.

### DICOM Server




Store pure and raw DICOM images to a DICOM device.

### Worklist/RemoteArch-RemoteHD




Search in a DICOM Modality Worklist, the patient found is copied into a remote database. The patient information and the examination results are stored to a remote database. Images are stored to an image network volume.

### RemoteArch-RemoteHD/eVue



A remote database is used for patient archiving. Images are stored to a network image volume and an MPEG exam is created in the configured destination.

### Worklist - DICOM Server




Search in the DICOM Modality Worklist. Images are stored to a DICOM Server.

# Dataflows


## Selections

### DICOM Query Retrieve - DICOM Server




Retrieve images and Private SR from a DICOM server. Images are stored to a DICOM server. Some of the measurements are stored if DICOM SR is turned on

### DICOM Print



Prints images to a DICOM printer via DICOM spooler.

### Worklist/DICOM Query Retrieve -DICOM Server




Search in a DICOM Modality Worklist, retrieve images and Private SR from a DICOM server. Images are stored to a DICOM server. Some of the measurements are stored if DICOM SR is turned on



# Dataflows


## Selections

### DICOM CD/DVD read




Read DICOM Media from the CD/DVD-drive. Read only dataflow, no data can be stored. Purpose is to open patient and view exam images

### DICOM USB Hardisk/Memstick Read




Read DICOM data from an USB hard disk or memory stick. Read only dataflow, no data can be stored. Purpose is to open patient and view exam images

### Misc. Export




### DICOM CD/DVD



Read and write to DICOM Media from the CD/DVD-drive.. Purpose is for transfer (import/export)

### DICOM USB Hardisk/Memstick




Read and write DICOM data to an USB hard disk or memory stick. Purpose is for transfer (import/export)

### No Archive




Enables to perform an examination without storing the data to any archive.

### CD/DVD



Read and write to database located on a CD/DVD-drive. This destination should be used when the transferred data is intended to be imported in a Vivid E9 or an EchoPAC Software Only with software version 113 or earlier

### USB Hardisk/Memstick




Read and write to database located on an USB hard disk or memory stick. This destination should be used when the transferred data is intended to be imported in a Vivid E9 or an EchoPAC Software Only with software version 113 or earlier



# Dataflows


## Hide/Unhide/Add dataflows

Select **Dataflow** on the **Connectivity** tab to work with dataflows



### Predefined dataflows

The default dataflows are found in the pulldown list



Per default, only a few dataflows are listed in the Search/Create Patient Screen. Unselect **Hidden** in the Connectivity tab to display more Dataflows

**Note!** For some dataflows the Hidden field is «greyed out», meaning they cannot be unhidden

**Note!** The Local Archive Int.HD - Shared do not exist on the Scanner

# DICOM Activities

## Vivid E80/E90/E95

There are six local real-world activities that occur in Vivid - **Exam Save, Echo, Worklist Query, Image Store/End Exam, Image Print** and **Exam Query/Retrieve**

**Exam save** initiates a connection with the DICOM SCP and transmits images and results to the DICOM SCP. If Storage Commitment is configured a commitment request will be sent for the images and results.

**Echo** initiates a connection with the DICOM SCP, posts a Verification request and closes the connection. It also responds to incoming Verification requests (for service use).

**Worklist Query** initiates a connection with the DICOM SCP, performs a query and retrieves the matching entries to the product.

**Image Store/End exam:** If Modality Performed Procedure Step is configured N-CREATE and N-SET messages will be sent for the exam.

**Image Print** will send images to a DICOM Print SCP.

**Exam Query/Retrieve** initiates a connection with the DICOM SCP, performs a query and retrieves selected examination.

# Compatibility

## Vivid E80/E90/E95 & EchoPAC SWO

When a new scanner (E95, E90, E80, S60, S70) arrives at a site:

If site has ImageVault + RDCM: Must upgrade **RDCM**. **RDCM** requires 64 bit OS with DirectX11 support.

This implies that **ImageVault must be IV 5.0.1** (Win Server 2008 R2)

If site has EchoPAC SWO:

Must upgrade to EchoPAC SWO 201.x.x. However there are no Dicom Share on EchoPAC SWO 201.x.x

If DICOM Share is needed, use EchoPAC SWO 113.x.x

A DirectX 11 compliant GPU is required for interactive review of 3D images and also recommended for reviewing 2D images. (Most graphics cards newer than 3 years old computer have this)

EchoPAC SWO, EchoPAC Plugin, RDCM, Vivid E95, E90, E80:

Software version: 201.x.x

Raw data version: 9001

Will not accept RAW data from:

Vivid 7 BT01, BT02, BT03.


Vivid 3

Vivid 4


# Worklist

## Configure


The DICOM Worklist is an **Input** device in the Workflow/Dataflow, and the Vivid must be registered as a user on the Worklist Server.



From dataflow name pulldown list, select the appropriate **Worklist** dataflow



Select **DicomWorklist** from the Dataflow Overview window



Select **Properties**

**Required fields**

The **Worklist** Properties window will appear. The required fields ★ must be filled out. Other fields ♦ may need to be adjusted

Enter correct **AE title**, as it is case sensitive and **Port No.** **IP Addr** must be selected from the pulldown menu and **Name** is optional, but recommended

If **IP Addr** is not in pulldown menu, go to **TCPIP** tab and select **Modify** to change or **Add**. A **Server Config** window will appear

Select **Add** to insert a new, select **Modify** to correct **IP Addr** of existing DICOM Device. Insert a descriptive **Name** and **IP Address**

**Check** will ping the server

# Worklist

## Search Criterias

The DICOM Worklist is an **Input** device in the Workflow/Dataflow, and the Vivid must be registered as a user on the Worklist Server. The modality for MWL in the past **must** be **US** to start an exam




The 'DicomWorklist - (SCBasicWorklistMan)' dialog box contains the following fields and buttons:

- IP-address: ( HL7 ) 10.0.0.7 (dropdown menu) and Check button
- Name: DicomWorklist01 (text field)
- AE Title: MERGE\_WORK\_SCP (text field)
- Port No: 107 (text field)
- Max. Result: 500 (text field)
- Search Criterias button
- Retry section with:
  - Max #: 0 (text field)
  - Interval: 1 (text field) sec.
  - Timeout: 30 (text field) sec.
- OK and Cancel buttons



Select **Search Criterias**



The 'Search Criterias' dialog box contains the following fields and buttons:

- Select Tag: 00400001 Scheduled Station AE Title (dropdown menu)
- Value: XXXXXXXX (text field) and Don't Use checkbox
- Update list and Remove buttons
- Table with columns: Name, Value, Don't Use
- Table content: 00080060 Modality US
- OK and Cancel buttons

**Search Criterias** can be used to narrow the search in a Worklist

A Search Criteria properties window appears. In the Select Tag pulldown menu select the search criteria. Enter a value in **Value** if necessary.

Select **Update list** to add search criteria. To remove, select the search criteria to delete from the **Name** listbox then select **Remove**

Select **OK** to close window

# Worklist

## Connect

Per default, Worklist dataflows are hidden, Uncheck the «Hidden» checkbox in the dataflow page to include the dataflow

- ☐ Default
- ☐ Direct Store
- ☒ Hidden



The dialog box contains the following fields and controls:

- IP-address: (DICOMSERVER) 10.0.0.5 [Check]
- Name: DicomWorklist01
- AE Title: MERGE\_WORK\_SCP
- Port No: 107
- Max. Result: 500
- Search Criteria: [button]
- Retry: [button]
- Max #: 0
- Interval: 1 sec.
- Timeout: 30 sec.
- OK [button]
- Cancel [button]


In the **Worklist** Properties window select **Check** to test the Worklist service



Worklist service is available at selected server




Server cannot provide Worklist service



The dataflow overview window shows a tree structure with 'DicomWorklist' selected under 'Inputs'. To the right, there is a 'Repeats:' field with the value '3' and a 'Check' button. Below these is a 'Properties' button.


From the dataflow overview window, select **Worklist**



The dialog box shows 'Repeats:' with the value '3' and a 'Check' button.


Select **Check** button

DICOM ping passed  
(3 attempts):



The dialog box shows 'Repeats:' with the value '3', a 'Check' button, and a green checkmark icon.

DICOM ping failed  
(3 attempts):



The dialog box shows 'Repeats:' with the value '3', a 'Check' button, and a red X icon.



The dialog box displays the following information:

- Server Check (Ping): 3 of 3 OK
- DICOM Ping: 0 of 3 OK
- Ok [button]

Check status window



### Server Check:

**X of X OK** – Server IP Address is correct

**0 of X OK** – Server IP Address is incorrect

### DICOM Ping:

**X of X OK** – Server **AE title** and/or **Port No** is correct

**0 of X OK** – Server **AE title** and/or **Port No** is incorrect

# DICOM Storage

The DICOM Storage is an **Output** device in the Workflow/Dataflow.

## Storage Device Configuration Scenarios

The VividE95/85/80 and EchoPAC may be connected to different storage device scenarios in a Hospital network, which will require different configurations, the four possible scenarios are:

**Scenario A:** Non Vivid/EchoPAC Dicom review station

**Scenario B:** Dicom review solution with EchoPAC Plugin

**Scenario C:** Mix of EchoPAC and Dicom review station

**Scenario D:** EchoPAC review station only


Storage Device	Dicom MF	Private SR data	RAW data
Scenario A	✓	✗	✗
Scenario B	✗	✓	✓
Scenario C	✓	✓	✓
Scenario D	✗	✓	✓

Storage Device Configuration Scenarios


**Note:** In a DICOM workflow, enable Private SR data to store all measurements and other exam data in the SR. This information is retrieved when the exam is loaded on the scanner, EchoPAC SWO and EchoPAC Plugin. In a DICOM workflow, this is the only way we can retrieve more than just images and the most basic exam information.

# DICOM Storage


## Configure



From dataflow name pulldown list, select the appropriate **DICOM** dataflow




Select **DicomStorage** from the Dataflow Overview window




Select **Properties**


Required fields




The **DicomStorage** Properties window will appear. The required fields ★ must be filled out. Other fields ◆ might need to be adjusted



Enter correct **AE title**, as it is case sensitive and **Port No.** **IP Addr** must be selected from the pulldown menu and **Name** is optional, but recommended



If **IP Addr** is not in pulldown menu, go to **TCPIP** tab and select **Modify** to change or **Add**. A **Server Config** window will appear



Select **Add** to insert a new, select **Modify** to correct **IP Addr** of existing DICOM Device  
Insert a descriptive **Name** and **IP Address**

**Check** will ping the server



# DICOM Storage Settings

If Frame rate, Compression and Quality settings are not set correctly, this may lead to long transfer time, huge file size and images not suited for diagnostics

**Note!**  
Image settings only apply to Dicom SF and DICOM MF.

Verify this default ★ Image Settings:  
**Only Black/White** - unchecked

Adjust these Image Settings if needed:

Max Framerate  
Compression  
Quality %  
Allow Raw Data  
Raw Compr.  
Allow Multiframe

25/Full  
JPEG  
95  
Adjust according to Scenario A,B,C or D

★ ☐ Reopen pr. Image

Verify this default ★ Image Setting:  
**Reopen pr. Image** – unchecked

If Reopen per image is enabled, VIVID will create a new connection (association) for each image. This may be useful for DICOM servers that do not accept different image types in the same association.

Verify these default ★ Dicom SR Settings:

**No Images** – set if SR data has to be sent without images (different server for Images)  
**Signed Doppler Velocities** – unchecked  
**Use older SR version** – unchecked

Adjust these DICOM SR Settings if needed:

**Allow SR** - set if DICOM supported M&As are to be sent to a DICOM SR Server  
**Allow SR Private Data** – check to allow Private tags

Adjust according to Scenario A,B,C or D

Timeout 40 sec.

Some DICOM Servers may require a higher Timeout, than default

# DICOM Storage

## Storage Commitment

Dicom Storage Commitment is a way for an Image Storage Service to inform your scanner that it has successfully stored your images without corruption.

☒ Storage commitment

Select **Storage commitment**

The **Storage Commitment** properties window will appear. The required fields ★ must be filled out. Other fields ⚡ might need to be adjusted

Enter correct **AE title**, as it is case sensitive and **Port No.** **IP Addr** must be selected from the pulldown menu and **Name** is optional, but recommended

If **IP Addr** is not in pulldown menu, go to **TCPIP tab** and select **Modify** to change or **Add**. A **Server Config** window will appear


**Check** will ping the server

Select **Add** to insert a new, select **Modify** to correct **IP Addr** of existing DICOM Device  
Insert a descriptive **Name** and **IP Address**

# DICOM Storage

## MPPS


A Dicom Performed Procedure, also known as Modality Performed Procedure Step (MPPS) is a way to inform your billing system that an ultrasound procedure has been completed and can be billed



The 'DicomStorage - (SCStorage)' window shows various settings. The 'Storage commitment' checkbox is highlighted with a red box, and the 'MPPS' checkbox is also highlighted. The 'Image settings' section includes 'Max Framerate' (Full), 'Compression' (Jpeg), and 'Quality %' (95). The 'Dicom SR Settings' section includes 'Allow SR', 'Allow SR Private Data', 'No Images', 'Signed Doppler Velocities', and 'Reopen pr. Image'. The 'Retry' section includes 'Max #' (1), 'Interval' (120 sec.), and 'Timeout' (40 sec.).



Select MPPS




The 'MPPS properties' window shows fields for 'IP-address' (DICOMSERVER) 10.0.0.5, 'Name' (DicomMpps01), 'AE Title', and 'Port No'. The 'Retry' section includes 'Max #' (1), 'Interval' (120 sec.), and 'Timeout' (40 sec.). The 'AE Title' and 'Port No' fields are marked with red stars, indicating they are required. The 'Timeout' field is marked with an orange diamond, indicating it might need to be adjusted.

The **MPPS** properties window will appear. The required fields ★ must be filled out. Other fields ♦ might need to be adjusted




The 'MPPS properties' window shows fields for 'IP-address' (DICOMSERVER) 10.0.0.5, 'Name' (DicomMpps01), 'AE Title', and 'Port No'. The 'Retry' section includes 'Max #' (1), 'Interval' (120 sec.), and 'Timeout' (40 sec.). The 'AE Title' and 'Port No' fields are marked with red stars, indicating they are required. The 'Timeout' field is marked with an orange diamond, indicating it might need to be adjusted.

Enter correct **AE title**, as it is case sensitive and **Port No.** **IP Addr** must be selected from the pulldown menu and **Name** is optional, but recommended



The 'MPPS properties' window shows fields for 'IP-address' (DICOMSERVER) 10.0.0.5, 'Name' (HCE-HNYKYW1) 127.0.0.1, 'AE Title' (ECHOPAC-000000) 10.0.0.4, and 'Port No' (DICOMSERVER) 10.0.0.5. The 'Retry' section includes 'Max #' (1), 'Interval' (120 sec.), and 'Timeout' (40 sec.).

If **IP Addr** is not in pulldown menu, go to **TCPIP** tab and select **Modify** to change or **Add**. A **Server Config** window will appear



The 'CONNECTIVITY' window shows the 'Tcpip' tab. The 'Server Config' section lists servers: (HCE-HNYKYW1) 127.0.0.1, (ECHOPAC-000000) 10.0.0.4, (DICOMSERVER) 10.0.0.5, (HL7) 10.0.0.7, (NAS) 192.168.1.60, and (GW) 192.168.1.1. The 'Add' button is highlighted. The 'Server Config' window shows fields for 'Server Name' (New\_Name) and 'IP-Address' (0.0.0.0).

Select **Add** to insert a new, select **Modify** to correct **IP Addr** of existing DICOM Device. Insert a descriptive **Name** and **IP Address**


**Check** will ping the server

# DICOM Storage

## Connect

Per default, DICOM Storage dataflows are hidden, Uncheck the «Hidden» checkbox in the dataflow page to include the dataflow

- ☐ Default
- ☐ Direct Store
- ☒ Hidden



The 'DicomStorage - (SCStorage)' window contains the following fields and options:


- IP-address: (DICOMSERVER) 10.0.0.5 (with a 'Check' button)
- Name: DicomStorage01
- AE Title: DICOMSTORAGESCP (with a 'Storage commitment' checkbox)
- Port No: 105 (with an 'MPPS' checkbox)
- Image settings:**
  - Max Framerate: Full
  - Compression: Jpeg
  - Quality %: 95
  - Allow Raw Data (checkbox)
  - Raw Compr. (checkbox)
  - Allow Multiframe (checkbox)
  - Only black/white (checkbox)
- Dicom SR Settings:**
  - Allow SR (checkbox)
  - Allow SR Private Data (checkbox)
  - No Images (checkbox)
  - Signed Doppler Velocities (checkbox)
  - Reopen pr. Image (checkbox)
- Retry:**
  - Max #: 1
  - Interval: 120 sec.
  - Timeout: 40 sec.
- Buttons: OK, Cancel



DICOM Storage service is available at selected server




Selected server cannot provide DICOM Storage service



The 'Dataflow overview' window shows a tree structure with 'DicomStorage' selected under 'Outputs'. To the right, there is a 'Repeats:' field set to 3, and 'Check' and 'Properties' buttons.

From the dataflow overview window, select **DicomStorage**




The 'Repeats:' dialog box shows the value 3 and a 'Check' button.

Select **Check** button


In the **DicomStorage** Properties window select **Check** to test the DICOM Storage service

DICOM ping passed (3 attempts):




The 'Repeats:' dialog box shows the value 3, a 'Check' button, and a green checkmark icon.

DICOM ping failed (3 attempts):



The 'Repeats:' dialog box shows the value 3, a 'Check' button, and a red X icon.



The 'Information - Check status' window displays the following information:

- Server Check (Ping): 3 of 3 OK
- DICOM Ping: 0 of 3 OK
- Buttons: Ok

Check status window

**Server Check:**

**X of X OK** – Server IP Address is correct

**0 of X OK** – Server IP Address is incorrect

**DICOM Ping:**


**X of X OK** – Server **AE title** and **Port No** is correct

**0 of X OK** – Server **AE title** and **Port No** is incorrect

# DICOM CD/DVD

## Configure

The **DICOM CD/DVD read** is an input dataflow and is to open patient and view exam images, while the **DICOM CD/DVD** is both input and output dataflow



From dataflow name pulldown list, select the appropriate **DICOM CD/DVD** dataflow

**DICOM CD/DVD read** dataflow:



Select **Dicom CD/DVD** from the Dataflow Overview window

**DICOM CD/DVD** dataflow:




Select **Dicom CD/DVD** from the Dataflow Overview window

Properties

Select **Properties**



The **DICOM Media** Properties window for the **DICOM CD/DVD read** dataflow will appear.



The **DICOM Media** Properties window for the **DICOM CD/DVD** dataflow will appear.

Configure - DICOM CD/DVD

# DICOM CD/DVD

## Settings

If Frame rate, Compression and Quality settings are not set correctly, this may lead to long transfer time, huge file size and images not suited for diagnostics

DicomMedia - (SCMediaStorage)

IP-address ( HCE-HNYKYW1 ) 127.0.0.1 Check

Name DICOM CD/DVD

Location \\127.0.0.1\CDRW

Type R/W

**Image settings**

Max Framerate Full ☐ Allow Raw Data

Compression Jpeg ☐ Raw Compr.

Quality % 95 ☒ Allow Multiframe ☐ Only black/white

**Dicom SR Settings**

☐ Allow SR ☐ Use older SR version

☐ Allow SR Private Data

☐ No Images

☐ Signed Doppler Velocities

OK Cancel

DICOM CD/DVD dataflow properties

**Image settings**

Max Framerate Full ☐ Allow Raw Data

Compression Jpeg ☐ Raw Compr.

Quality % 95 ☒ Allow Multiframe ☐ Only black/white

Verify these default ★ Image Settings:

**Allow Raw Data** – unchecked  
**RAW Compr** – greyed out  
**Allow Multiframe** – checked  
**Only Black/White** – unchecked

Adjust these ⚡ Image Settings if needed:

**Max Framerate** → 25/Full  
**Compression** → JPEG  
**Quality %** → 95

**Dicom SR Settings**

☐ Allow SR ☐ Use older SR version

☐ Allow SR Private Data

☐ No Images

☐ Signed Doppler Velocities


Verify these default ★ Dicom SR Settings:

**Allow SR** – set if DICOM supported M&As are to be sent to a DICOM SR Server  
**Allow SR Private Data** – greyed out  
**No Images** – set if SR data has to be sent without images  
**Signed Doppler Velocities** – greyed out  
**Use older SR version** – greyed out

# DICOM USB

## Configure

The **DICOM USB read** is an input dataflow, and is to open patient and view exam images, while the **DICOM USB** is both input and output dataflow



From dataflow name pulldown list, select the appropriate **DICOM USB** dataflow

**DICOM USB read** dataflow:




Select **DICOM USB** from the Dataflow Overview window

**DICOM USB** dataflow:




Select **DICOM USB** from the Dataflow Overview window

Select **Properties**



The **DICOM Media** Properties window for the **DICOM USB read** dataflow will appear.



The **DICOM Media** Properties window for the **DICOM USB** dataflow will appear.

# DICOM USB

## Settings

If Frame rate, Compression and Quality settings are not set correctly, this may lead to long transfer time, huge file size and images not suited for diagnostics

DicomMedia - (SCMediaStorage)

IP-address (AURORA-000132 ) 127.0.0.1 Check

Name DICOM USB Harddisk/Memstick

Location \\127.0.0.1\MEMSTICK

Type R/W

**Image settings**

Max Framerate Full ☐ Allow Raw Data

Compression Jpeg ☐ Raw Compr.

Quality % 95 ☒ Allow Multiframe ☐ Only black/white

**Dicom SR Settings**

☐ Allow SR ☐ Use older SR version

☐ Allow SR Private Data

☐ No Images

☐ Signed Doppler Velocities

**Dicom Viewer**

☐ Add Dicom Viewer

OK Cancel

DICOM USB dataflow properties

**Image settings**

Max Framerate Full

Compression Jpeg

Quality % 95

Verify these default ★ Image Settings:

**Allow Raw Data** – unchecked  
**RAW Compr** – greyed out  
**Allow Multiframe** – checked  
**Only Black/White** - unchecked

Adjust these Image Settings if needed:

**Max Framerate** → 25/Full  
**Compression** → JPEG  
**Quality %** → 95

**Dicom SR Settings**

☒ Allow SR ☐ Use older SR version

☒ Allow SR Private Data

☒ No Images

☒ Signed Doppler Velocities

Verify these default ★ Dicom SR Settings:


**Allow SR** – set if DICOM supported M&As are to be sent to a DICOM SR Server  
**Allow SR Private Data** – greyed out  
**No Images** – set if SR data has to be sent without images  
**Signed Doppler Velocities** – greyed out  
**Use older SR version** – greyed out




# Query Retrieve

## Configure

The DICOM Query Retrieve is an **Input** device in the Workflow/Dataflow, and the Vivid must be registered as a user on the Q/R DICOM Server. Dicom QR configuration on EchoPAC might also require Firewall modifications




From dataflow name pulldown list, select the **DICOM Query retrieve** dataflow



The **Query Retrieve** Properties window will appear. The required fields ★ must be filled out. Other fields ◆ might need to be adjusted

Enter correct **AE title**, as it is case sensitive and **Port No.** **IP Addr** must be selected from the pulldown menu and **Name** is optional, but recommended




Select **Add** to insert a new, select **Modify** to correct **IP Addr** of existing DICOM Device  
Insert a descriptive **Name** and **IP Address**

**Check** will ping the server

# Query Retrieve

## Search Criterias

The DICOM Query Retrieve is an **Input** device in the Workflow/Dataflow, and the Vivid must be registered as a user on the Q/R DICOM Server.




The 'DicomQR - (SCQueryRetrieve)' dialog box contains the following fields and controls:

- IP-address: ( DICOMSERVER ) 10.0.0.5 (with a dropdown arrow) and a 'Check' button.
- Name: DICOM Query retrieve
- AE Title: AE\_DICOMSERVER
- Port No: 110
- Max. Result: 100
- A 'Search Criterias' button.
- A 'Retry' section with:
  - Max #: 0
  - Interval: 1 sec.
  - Timeout: 20 sec.
- 'OK' and 'Cancel' buttons at the bottom.



**Select Search Criterias**



The 'Search Criterias' dialog box contains the following fields and controls:

- 'Select Tag' dropdown menu: 00080030 Study Time
- 'Value' text field and a 'Don't Use' checkbox.
- 'Update list' and 'Remove' buttons.
- A listbox with the following items:

Name	Value	Don't Use
00080061 Modalities in Study	US	
00080060 Modality		
- 'OK' and 'Cancel' buttons at the bottom.

**Search Criterias** can be used to narrow the search in a Query Retrieve

A Search Criteria properties window appears. In the Select Tag pulldown menu select the search criteria. Enter a value in **Value** if necessary.

Select **Update list** to add search criteria. To remove, select the search criteria to delete from the **Name** listbox then select **Remove**

Select **OK** to close window

# Query Retrieve

## Connect

Per default, QueryRetrieve dataflows are hidden, Uncheck the «Hidden» checkbox in the dataflow page to include the dataflow

- ☐ Default
- ☐ Direct Store
- ☒ Hidden



The dialog box contains the following fields and controls:

- IP-address: ( DICOMSERVER ) 10.0.0.5 [Check]
- Name: DICOM Query retrieve
- AE Title: AE\_DICOMSERVER
- Port No: 110
- Max. Result: 100
- Search Criteria: [button]
- Retry section: Max # 0, Interval 1 sec, Timeout 20 sec.
- OK and Cancel buttons.


In the **DicomQR** Properties window select **Check** to test the QueryRetrieve service



DicomQR service is available at selected server




Selected server cannot provide DicomQR service



The window shows a dataflow diagram with 'DicomQR' in the center. To the right, there is a 'Repeats' field set to 3, a 'Check' button, and a 'Properties' button.


From the dataflow overview window, select **DicomQR**



The dialog box shows 'Repeats: 3' and a 'Check' button.


Select **Check** button

DICOM ping passed (3 attempts):




The dialog box shows 'Repeats: 3' and a 'Check' button with a green checkmark icon next to it.

DICOM ping failed (3 attempts):



The dialog box shows 'Repeats: 3' and a 'Check' button with a red X icon next to it.



The dialog box displays the following information:

- Server Check (Ping): 3 of 3 OK
- DICOM Ping: 0 of 3 OK
- Ok button

Check status window

**Server Check:**

**3 of 3 OK** – Server IP Address is correct

**0 of 3 OK** – Server IP Address is incorrect

**DICOM Ping:**


**3 of 3 OK** – Server **AE title** and **Port No** is correct

**0 of 3 OK** – Server **AE title** and **Port No** is incorrect

# DICOM Spooler

## Check Status

To check status of the DICOM Spooler, select **F4** or **Alt** + **S** on the alphanumeric keyboard or select **Spooler** in the Utility menu on the Touchpanel.



Touch Screen

DICOM spooler displays the current DICOM output jobs, and is used for checking the current job's status when a job is saved or displays an error



DICOM Job Spooler window

To remove the DICOM Spooler Job window select


Hide

# DICOM Printer

## Modify Outputs


Select **Additional Outputs** tab to get to the Additional Outputs page for configuring the **P1** and **P2** buttons on the Operator Panel.

Per default **P1** is set to Printer,  
**P2** is set to store to clipboard




The **P1** and **P2** position on the Operator panel

Select the appropriate button from the pulldown menu



Additional Outputs window


### Add:



In **Available outputs** select the device you want to use as output




Select > to add the output device



Verify that the device you selected, appears in the **Selected devices** window

**Note !** Using a Print button for image transfer to a Dicom service is a quick way to accomplish "direct store" to that service without having to set up a dataflow


### Remove:



In **Selected devices** window select the device you want to remove




Select < to remove the output device




Verify that the device you selected, disappears in the **Selected devices** window

# DICOM Printer

## Configure




Select **Dicom Print** in the Selected Output window




Select **Properties...**


Required fields




The **DICOM Print** Properties window will appear. The required fields ★ must be filled out. Other fields ♦ might need to be adjusted



Enter correct **AE Title**, as it is case sensitive and **Port No.** **IP Addr** must be selected from the pulldown menu and **Name** is optional, but recommended



If **IP Addr** is not in pulldown menu, go to **TCPIP** tab and select **Modify** to change or **Add**. A **Server Config** window will appear




Select **Add** to insert a new, select **Modify** to correct **IP Addr** of existing DICOM Device. Insert a descriptive **Name** and **IP Address**

**Check** will ping the server

# DICOM Printer

## Connect



The DicomPrint - (SCPrintManagement) dialog box is shown. It has tabs for Densities, Settings, and Job. The IP-address field is set to (HCE-HNYKYW1) 127.0.0.1, and the Name field is set to DicomPrint02. The AE Title and Port No fields are empty. The Interval is set to 120 sec, and the Timeout is set to 40 sec. The Max # field is set to 1. The Densities tab is selected, showing Min and Max density fields set to 0, and Border and Empty image dropdowns set to BLACK and WHITE respectively. OK and Cancel buttons are at the bottom.




DicomPrint service is available  
at selected device



Selected device cannot  
provide DicomPrint service


In the **DicomPrint** Properties  
window select **Check** to test the  
DicomPrint service

If DicomPrint service failed or printout is not satisfying, verify the service configurations in the Densities, Settings and Job tabs:



The Densities tab of the DicomPrint dialog is shown. It contains fields for Min and Max density, both set to 0. There are dropdown menus for Border (set to BLACK) and Empty image (set to WHITE). OK and Cancel buttons are at the bottom.

Densities



The Settings tab of the DicomPrint dialog is shown. It contains fields for Format (STANDARD1.1), Film Orientation (LANDSCAPE), Magnification Type (REPLICATE), Film Size (8INX10IN), Smoothing Type, and Trim (YES). There is a Configuration Information section at the bottom. OK and Cancel buttons are at the bottom.

Settings



The Job tab of the DicomPrint dialog is shown. It contains fields for Film session label (1), Film destination (MAGAZINE), Medium (PAPER), Color (GREY), Priority (HIGH), and Copies (1). OK and Cancel buttons are at the bottom.


Job

# Network Printer


## Add New

**Note!** Ensure the Network Printer is powered on and connected to the site network before adding the Printer

Select  on the **Connectivity** tab to add Network or Local Printers

Select  on the **Service menu** to add a Printer

A Standard Windows Printer installation program will start,



Select this option to install a Local Printer

Select this option to install a Network Printer

Window Printer Installation

Click  to continue

Follow the instructions to add a network printer successfully




# Network Printer

## Configure IP Address

**Note!** Use this method if Printer is already installed, but Printer IP Address needs to be reconfigured

Select **Service** on the **Connectivity** tab to configure Network Printers

In the **Network Printer** section enter the correct IP address of the Network Printer



**Network Printer** section

When the new IP address is set, the **Set IP-address** button is enabled. Select this button to store the new IP address to the Network Printer

Select the correct name for the Network Printer, if no name is correct select the first Generic Name available


# Network Printer

## Rename


**Note!** Use this method if Printer is already installed, but Printer Name needs to be reconfigured

Select **Additional Outputs** the **Connectivity** tab to rename Network Printers


Select **Queue** in the **Printer Setup** section for the Printer Queue window to appear



Printer Setup section



Select **Printer** in the **Printer Queue** window to get to **Printer Properties**



Printer Pull down menu


Printer Queue window

Select **Printing Preferences** to adjust advanced printer settings



Select **Properties** to change the Printer name


Enter the new Network Printer name here. Select **Apply** and click **OK** to store the new Printer Name




Printer Properties window

# Remote Archive


## Configure




From dataflow name pulldown list, select the appropriate **Remote Archive** dataflow




Select **Cardiology\_3\_1\_4\_ArchiveService** from the Dataflow Overview window



Select **Properties**




The **Cardiology\_3\_1\_4\_ArchiveService** Properties window will appear. If DNS is not set up on the network, select the **Update 'hosts' file** button.  
**Note!** On EchoPAC SWO, Windows administrative rights are needed to do this



Enter correct **IP-address** of the Remote Archive server. The other fields will be populated automatically.



If **IP Addr** is not in pulldown menu, go to **TCPIP** tab and select **Modify** to change or **Add**. A **Server Config** window will appear




Select **Add** to insert a new, select **Modify** to correct **IP Addr** of Remote server

Insert correct **HostName** and **IP Address**

**Check** will ping the server

# Remote Archive

## Connect




In the **Cardiology\_3\_1\_4\_ArchiveService** Properties window select **Check** to ping the remote server




Remote server IP Address is correct



Remote server IP Address is incorrect




From the dataflow overview window, select **Cardiology\_3\_1\_4\_ArchiveService**




Select **Check** button


Server Database and Image archive are correct (3 attempts):



Server Database and Image archive are incorrect (3 attempts):



### Check status window



### Server Check:

**X of X OK** – Server IP Address is correct  
**0 of X OK** – Service IP Address is incorrect

### Database Check:

**X of X OK** – Remote DataBase is found  
**0 of X OK** – Remote DataBase is not found

### File Destination Check:

**X of X OK** – Image Shared folder is found  
**0 of X OK** – Image Shared folder is not found

# Backup & Restore

## Backup

Select **Backup** tab from the Admin menu

Archive to backup	Result	Last successful backup
<input type="checkbox"/> Patient Archive		No record
<input type="checkbox"/> System Configuration		No record

Destination Device: CD/DVD Writable (G:)

Remote Path:

Start backup

Backup window

Archive to backup	Result	Last successful backup
<input checked="" type="checkbox"/> Patient Archive		No record
<input checked="" type="checkbox"/> System Configuration		No record

Destination Device: USB HD/Memstick (I:)

Remote Path:

Start backup

- The media used for backup can not be used for storing images
- When Remote Path is selected, the system will use the default user/password or the configurable user/password from the **Tcpip** tab

Start backup

Select  
Start backup

Information

Backup is completed

Ok


Backup completed

Select **Patient Archive** and **System Configuration** and the correct **Destination Device** and/or **Remote Path**


# Backup & Restore

## Restore

Select **Restore** tab from the Admin menu



Restore window




Restore Now

Select  
Restore Now





To start **Restore** select **Admin** from the Config menu

From the **Restore** window there is a possibility to restore to Factory defaults configuration. Select the appropriate button and press **Restore Now**



Press **Ok** to confirm replacing settings to **Factory defaults**

The Restore functionality requires acknowledgement from the user to restore previous backup data, thus the similar pop-up windows will appear



The system will shut down after Restore is completed

**Note!** Backup & Restore may not work across all CVUS products and BTs..

Select **Patient Archive** and/or **System Configuration** and the correct **Destination Device** and **Remote Path**

# Disk Management Configuration

To get to **Disk Management** select the Connectivity menu

Select **Disk Management** tab from the Connectivity menu

## Predefined Reminder Interval

Set the reminder interval for Disk Management. The Recommended Interval is 1 week

## Manage files Older Than

Select the time for how old your image files should be before a backup is required. Recommended Interval is 1 week

## Operation

Select between 3 types of operations:  
Copy, Move and Delete  
Recommended operation is Move

The screenshot shows the 'CONNECTIVITY' menu with the 'Disk Management' tab selected. Red boxes highlight the following elements: the 'Reminder interval' dropdown set to '4 Weeks', the 'Manage files older than' dropdown set to 'Today (all files)', the 'Operation' radio buttons with 'Copy' selected, the 'Destination device' dropdown set to 'USB/HD', and the 'Run' button. Red lines connect these highlighted elements to their respective descriptions on the left side of the slide.

A legend showing the three operation options: 'Copy' (selected with a filled radio button), 'Move' (unselected with an empty radio button), and 'Delete' (unselected with an empty radio button).

**Copy** – Images and Reports will be copied to the Destination Device

**Move**– Images and Reports will be moved to the Destination Device

**Delete**– Images and Reports will be deleted.

## Run

Select «Run» to start Disk Management

## Destination Device

Select a removeable media or a network share folder

# Disk Management


## Running

Disk Management requires an acknowledgement from User that the correct settings have been selected and the calculated data size is displayed. Press Start to run Disk Management



Disk Management copies/moves the files found for the selected «Older Than» time window to the Destination device. There's a possibility to Abort the operation and resume later.

### Note!

After Disk Management is finished a Backup of the Patient Archive should be initiated



When Disk Management process completes, a summary status is displayed. There's a possibility to print an Exams list, Media list or a complete summary.






# EchoPAC SWO

## Local Archive Setup


**Note !** Per default the **Local Archive Int. HD - Shared** dataflow is hidden. To avoid confusion do not unhide this, unless the intention is to set up as an EchoServer




Released

# MPEGVue Export

## Configure





From dataflow name pulldown list, select the appropriate **Misc Export** dataflow



Select **MPEGVueService** from the Selected Outputs



Select **Properties**




Select the **Destination** from the drop down menu and/or enter Remote Path


Optionally Select to **Include Reports**

# MPEGVue Export


## Connect



In the **Selected devices** window select **MPEGvue**



Select **Properties** button



If **Remote path** was selected  
**Check** result:



Successfully logged into remote destination, write access confirmed



User and/or password incorrect or no write access to destination device




Go to «Archive» screen and select




to initiate MPEGVue file transfers

# XML Export


## Configure




From dataflow name pulldown list, select the appropriate **Misc Export** dataflow



Select **XmlExportService** from the Selected Outputs




Select **Properties**




Select the **Destination** from the drop down menu and/or enter Remote Path

# XML Export


## Connect



In the **Selected devices** window select **XmlExportService**



Select **Properties** button



If **Remote path** was selected  
**Check** result:



Successfully logged into remote destination, write access confirmed



User and/or password incorrect or no write access to destination device



Go to «Archive» screen and select




to initiate XmlExport file transfers

# ViewPoint

## Configure

### From ViewPoint:

Obtain the ViewPoint **IP-Address**, **AE-Title** and **Port No**,  
The **AE-Title** and **Port No.** is available on ViewPoint at the ViewPointAdm in the **Service Configuration** window



**Note!** Each ViewPoint has its own Unique AE-Title. The AE-Title displayed is for illustration only

### On Scanner:

Select the **Tcpip** to ensure the ViewPoint Server has been added to the Servers list in the **Server Config** section

If not present, select **Add** to add the ViewPoint server, then go to the **Dataflow** tab and select the appropriate dataflow

If present, go to the **Dataflow** tab and select the appropriate dataflow

If Worklist, DICOM Query/Retrieve was enabled when configuring the ViewPoint the appropriate dataflow is :

[Worklist/DICOM Query Retrieve -DICOM Server](#)

# ViewPoint Connect

To connect the Scanner to ViewPoint, the ViewPoint must be configured first. On ViewPoint, select the **ViewPointAdmin**



Log in with User name and Password

Log on to VPAdmin


User name:

Password:

OK Cancel

VPAdmin Log in

Select **Services Manager** in VPAdmin window




Select the **DICOM Server** in the services list

VPAdmin window, Service Manager tab

Select **Stop** to stop the **DICOM Server** service


Select **Configure...** to configure the **DICOM Server** client.

A Service configuration window will appear



Service Configuration window

Select **Add...** to add a new **DICOM Server** client. The **Modality Setting** window will appear




Modality Setting window

Select **OK** twice to close both the **Modality Setting** and **Service Configuration** window and **Save now** on the pop-up dialog

# ViewPoint

## Connect cont.

Select **ID Mapping** tab in **VPAdmin** window




VPAdmin window, ID Mapping tab

Check **Ext. ID 2** ☒ and select **Options...**



The **External ID 2** Properties window will appear




External ID 2 Properties window

From the pull down menu select **DICOM Patient ID**

Select **OK** when done

Select **File** in the **VPAdmin** menu bar



From the **File** pull down menu select **Save** to save the changes

Select **Services Manager** tab in **VPAdmin** window

Select the **DICOM Server** in the services list

Select **Start** to start the **DICOM Server** service



# Insite ExC

## Connect

### Pre-requisites & Preparation:

- To configure Insite ExC you need to know the rotating **GE service password**.
- Network setup complete, refer to «Network Settings» section
- Verify for VividE9 an **HTTPS (Secure Sockets Layer)** internet connection via port **443**
- If **proxy server** is used then ask Hospital IT for IP address, port number and proxy user and password.
- **System ID** (MUST-CARES-CAPS-Siebel)


The image displays three screenshots from the Insite ExC software interface, illustrating the steps to connect to the service. Red arrows and boxes highlight the specific elements mentioned in the numbered instructions.

- 1. Select Insite icon**: The first screenshot shows the top navigation bar with several icons. The 'Insite' icon, which features a telephone handset, is highlighted with a red box.
- 2. Select «Service Desktop»**: The second screenshot shows the 'SEARCH / CREATE PATIENT' screen. Within the 'InSite ExC' section, the 'Service Desktop' button is highlighted with a red box.
- 3. «Select User Level: GE Service»  
Enter rotating GE service password  
Select «Okay»**: The third screenshot shows the 'Service Login' dialog box. The 'Select User Level' dropdown menu is set to 'GE Service', and the 'Enter Password' field contains a series of dots. Both the dropdown and the password field are highlighted with a red box.

# Insite ExC

## Configure

4.  
Select «Configuration»



5.  
Select «Insite ExC Agent Configuration»

6.  
Complete the Configuration page, all items in **bold** are mandatory

The screenshot displays the 'Agent Configuration' form. Fields are highlighted with red boxes and labeled with red arrows. The form includes sections for general information, location, and server settings. The 'Enterprise Server' section is expanded, showing fields for URL and Tunnel. The 'Proxy Configuration' section is also visible at the bottom.

CRM No: System ID

Continent & Country

City & State(Prov)

Institution

**Enterprise Server:** Select **OTHER** and type in the URL addresses:

Enterprise Server: <https://198.169.188.10:443>


Enterprise Tunnel: <https://198.169.188.11:443>

**Proxy:** Only needed if customer uses Proxy Server

7.  
Select «Submit Changes»


# Insite ExC

## Configure Cont.



8.

«New agent configuration» screen will appear



9.

Close Service Desktop and Reboot the system. Repeat steps 1-3

10.

Select «Home» screen  
Service Connectivity Status will change to «Checked Out»

11.

All new connected systems to Insite should be reported to your OLC to initiate a «Device Grouping Request»

# Insite ExC

## RFS User Configuration

These Configuration steps are needed for the Auto Request For Service feature and must be repeated after each software update

1. Select Insite icon
2. Select «Request for Service»
3. Select «Users» and «Add User» in the «Permanent Users» section

The screenshot displays the Insite ExC software interface. The top menu bar includes 'Archive', 'Config', 'Playback', 'Help', and 'Exit'. Below the menu is a toolbar with buttons for 'New Patient', 'Open Exam', 'Export', 'Import', 'Move Exam', 'Delete', 'Print Patients', 'Clear Search', and 'Disk Manag'. The main window is titled 'SEARCH / CREATE PATIENT' and contains fields for 'Last Name', 'First Name', 'Category' (set to 'Cardiac'), 'Echolab', and 'Diagn. code'. A 'Patient ID' field is also present. A 'Request for Service' button is highlighted in the 'Service' section. The 'Contact GE - Service Browser' window is open, showing a 'Users' tab with a table of 'Permanent Users' and a table of 'Recent Users'. The 'Add User' button is highlighted in the 'Permanent Users' section.

# Insite ExC

## RFS User Configuration

Contact GE - Service Browser

Contact GE Queue Machine Queue Users

**Permanent Users**

Last	First	Phone	Ext.	Set D
GE	TEST	123456		

Add User

\* Last: GE  
\* First: TEST  
\* Phone: 123456  
Ext.:  
E-mail:

Add User

**Recent Users**

Last	First	Phone	Ext.
------	-------	-------	------

Make Selected Users Permanent Remove Selected Users

4. Complete all records marked with \*

5. Select «Add User»

6. Verify that there is a User marked in red in the «Permanent Users» section. If not, activate the desired user in the list and select «Set Default Machine Contact»

Contact GE - Service Browser

Contact GE Queue Machine Queue Users

**Permanent Users**

Last	First	Phone	Ext.	E-mail
GE	TEST	123456		

Add User Remove Selected Users Set Default Machine Contact

**Recent Users**

Last	First	Phone	Ext.	E-mail
------	-------	-------	------	--------

Make Selected Users Permanent Remove Selected Users

# Quick Aid

## Checklist



### Power

Verify that all connected devices are powered on.



### Connection

Verify that all connected devices have an Ethernet cable connected.



### Configuration

Verify that all connected devices have been configured correctly, with correct IP Address, Ports and AE titles.

# Quick Aid


## Troubleshooting Guide


### Network


#### Note!

Make sure to apply the GEHC Global Service Privacy and Security Standards to protect Sensitive Personal Information. (DOC1487129)

Before performing advanced troubleshooting, always assure that the following steps have been checked:



Regularly check the Network Status Icon at the lower left corner of the screen. When the Network Status Icon shows the following status :  The network connection is OK.

Symptoms	Possible Causes	Possible Solutions
<ul style="list-style-type: none"> <li>- Spooler full</li> <li>- Not able to export or import from a remote device</li> <li>- Not able to store images to a remote device</li> <li>- Timeout when connecting to a Remote Workflow</li> <li>- The Network Status icon shows the following status: </li> </ul>	<ul style="list-style-type: none"> <li>- Network cable disconnected, bad or loose</li> <li>- Hospital network down</li> <li>- Incorrect Network settings</li> <li>- Defective BEP</li> </ul>	<ul style="list-style-type: none"> <li>- Reseat Network cable</li> <li>- Replace Network cable</li> <li>- Verify TCP/IP settings in Connectivity page</li> <li>- Verify «Troubleshooting Remote Archive» and «Troubleshooting EchoPAC » sections</li> <li>- Contact Hospital IT for support</li> </ul>
<ul style="list-style-type: none"> <li>- Intermittent slowdown/lockup/freeze</li> <li>- System halts/pausing</li> <li>- Erratic movements of Cine loops</li> </ul>	<ul style="list-style-type: none"> <li>- Unstable Network</li> <li>- Slow network</li> <li>- Reading data from a Remote device on a slow network</li> </ul>	<ul style="list-style-type: none"> <li>- Contact IT for support and monitor Network</li> <li>- Contact IT to verify network speed</li> </ul>



# Quick Aid

## Troubleshooting Guide

### DICOM

#### Note!

Make sure to apply the GEHC Global Service Privacy and Security Standards to protect Sensitive Personal Information. (DOC1487129)


#### DICOM


☐ Detailed DICOM Log

#### Detailed DICOM Log

Always enable this for DICOM related problems, this will add detail DICOM information to the errorlogs. Found in the **Connectivity** menu in the **Tcpip** tab

Before performing advanced troubleshooting, always assure that the following steps have been checked:



Regularly check the Network Status Icon at the lower left corner of the screen. When the Network Status Icon shows the following status :  The network connection is OK.

Symptoms	Possible Causes	Possible Solutions
<ul style="list-style-type: none"> <li>- System images update freezes/halting</li> <li>- Not able to export to a DICOM Server</li> <li>- Not able to store loops in DICOM</li> <li>- Timeout when connecting to DICOM Device</li> <li>- Not able to transfer measurements to a DICOM Device</li> </ul>	<ul style="list-style-type: none"> <li>- DICOM Spooler full</li> <li>- Incorrect DICOM connectivity settings</li> <li>- Incorrect Workflow used</li> <li>- Slow Network speed</li> <li>- Incorrect DICOM SR configuration</li> </ul>	<ul style="list-style-type: none"> <li>- Empty DICOM Spooler.</li> <li>- Verify that the DICOM Device IP address, Port nr. and AE Title are correct</li> <li>- Change to Correct DICOM Workflow</li> <li>- Adjust timeout in DICOM connectivity settings</li> <li>- Adjust DICOM SR settings for the DICOM Device</li> </ul>
<ul style="list-style-type: none"> <li>- Not able to load DICOM Image</li> <li>- Images are not stored correctly</li> </ul>	<ul style="list-style-type: none"> <li>- DICOM dir may be corrupt</li> <li>- Incorrect DICOM connectivity settings</li> </ul>	<ul style="list-style-type: none"> <li>- Perform a repair DICOM DIR</li> <li>- Verify that DICOM Image settings are correct</li> </ul>




# Quick Aid


## Troubleshooting Guide

### Printer

Symptoms	Possible Causes	Possible Solutions
<ul style="list-style-type: none"> <li>- Unable to connect</li> <li>- Reports not printing</li> </ul>	<ul style="list-style-type: none"> <li>- Network cable disconnected or loose</li> <li>- Printer IP Settings are incorrect</li> <li>- Wrong Printer selected</li> <li>- Printer Driver Corrupt</li> <li>- Printer defective</li> </ul>	<ul style="list-style-type: none"> <li>- Reseat Network cable</li> <li>- Replace Network cable</li> <li>- Configure Printer with correct IP Address</li> <li>- Verify default printer is correct</li> <li>- Restore Printer Driver</li> <li>- Reload Software</li> <li>- Replace Printer</li> </ul>
<ul style="list-style-type: none"> <li>- Bad image printout</li> <li>- Report printout incorrect</li> </ul>	<ul style="list-style-type: none"> <li>- Low Ink</li> <li>- Printer Heads</li> <li>- Printer Layout Settings are incorrect</li> <li>- Corrupt report templates</li> <li>- Corrupt Printer Driver</li> <li>- Defective Printer</li> </ul>	<ul style="list-style-type: none"> <li>- Replace Ink Cartridges</li> <li>- Replace printer</li> <li>- Verify and adjust Printer Layout settings</li> <li>- Redo the Report template</li> <li>- Restore Printer Driver</li> <li>- Reload Software</li> <li>- Replace Printer</li> </ul>

Before performing advanced troubleshooting, always assure that the following steps have been checked:



Regularly check the Network Status Icon at the lower left corner of the screen. When the Network Status Icon shows the following status :  The network connection is OK.

# Quick Aid

## Troubleshooting Guide


### Note!


Make sure to apply the GEHC Global Service Privacy and Security Standards to protect Sensitive Personal Information. (DOC1487129)

### Archive

Symptoms	Possible Causes	Possible Solutions
<ul style="list-style-type: none"> <li>- Database corrupt</li> <li>- «No Archive»</li> <li>- Can not access internal Archive</li> <li>- Unable to log on</li> <li>- Not Archiving</li> </ul>	<ul style="list-style-type: none"> <li>- Database corrupt due to improper shutdown or sudden Power Loss</li> <li>- Harddrive full</li> </ul>	<ul style="list-style-type: none"> <li>- If no images stored, back up user presets, reload SW</li> <li>- Contact Service</li> <li>- Perform a Disk Management Move</li> </ul>
<ul style="list-style-type: none"> <li>- Not able to store to Remote Archive</li> <li>- Not able to store to DVD</li> </ul>	<ul style="list-style-type: none"> <li>- Network down</li> <li>- Incorrect Network settings</li> <li>- Incorrect CD/DVD disks used. DVD Drive defect</li> </ul>	<ul style="list-style-type: none"> <li>- Verify «Troubleshooting Remote Archive» and «Troubleshooting EchoPAC» sections</li> <li>- For Image archiving use CD/DVD -R/-RW disks</li> <li>- Replace DVD drive</li> </ul>

Before performing advanced troubleshooting, always assure that the following steps have been checked:




Regularly check the Network Status Icon at the lower left corner of the screen. When the Network Status Icon shows the following status :  The network connection is OK.

# Quick Aid

## Troubleshooting, Remote Archive

**Note!** EchoPAC sw version must be equal or higher than scanner sw version. ImageVault must be version 5.0.1.

**Database Check:** remote Sybase SQL Database of the EPPC/IV which listens for incoming Database traffic on TCP Port 2638.



**File Destination Check:** connection to the shared folder of the EPPC/IV which listens to TCP Port 445 for file Transfer (Image Files and Report Files).

## Troubleshooting – Remote Archive

### Database Check

#### Possible causes for failing:

- Remote Database not running.
- Too many DB connections - check “EchoShare” option on EPPC.
- Firewall closing inbound port 2638.
- Be sure you have compatible versions of EP or Image Vault software. EchoPAC must be at a version as high or higher than the version on the Vivid.

### File Destination Check

#### Possible causes for failing:

- File and Printer Sharing is disabled on the Remote Archive.
- Shared Folder not properly shared for users EchoClient and E1c2h3o4C5l6i7e8n9t, make sure that “Simple File Sharing” is turned off on the Remote Archive.
- Firewall on Remote Archive preventing File Sharing (make sure that Firewall allows file sharing and inbound port 445).

# Quick Aid

## Troubleshooting, EchoPAC

1.

Refer to «Readme» file on the installation media for minimum requirements like supported operating system, memory requirements, display settings, supported graphic boards. EchoPAC SWO requires **Win7 32 bits (or 64 bits). Win 8** is ok. Resolution of minimum **900 (in height)**. EchoPAC SWO requires **Direct3D Feature Level 11**. (Most graphics cards newer than 3 years old computer have this)

2.

The application must be installed with Administrator rights. Start EchoPAC SWO once with administrative rights after installation. If not, Local Archive will not work.

**Note:** On Windows 7 or newer: EchoPAC will not as default run with administrative rights even if user has administrative rights. To run EchoPAC with administrative rights: hold shift button, right-click on icon, select "Run As Administrator"

3.

Remote shared folder connectivity issues. Make sure that the share is not mapped on Windows level or disconnect all existing shares on Windows level and do a complete restart of Windows. Disable the "User Account Control" UAC before the installation of EPPC SWO

4.

If the PC is a member of Domain, the Domain Administrator must join the PC to a Workgroup before EPPC SWO installation

5.

Also refer to section Troubleshooting - Remote Archive:

# Quick Aid

## Repair DICOM Dir


This process is useful when copying multiple DICOM images to a read/writeable media, If the DicomDIR file is missing or corrupted, the DICOM Images are unreadable

Select **Tools** on the **Connectivity** tab to repair DicomDIR



Select the appropriate USB Memorystick or RW CD/DVD, from the Media pulldown menu. If media do not show, select **Refresh** button

If no DicomDIR was found the field DicomDIR is set to **No**



Select the **Repair DicomDIR** button, and verify a DicomDIR file is created

Refresh screen by selecting **Refresh**

Verify the field DicomDIR has changed and is set to **Yes**



# Quick Aid

**Note!** A corrupt archive may occur as a result of a improper or abnormal shut down.

## Local Archive Database Corrupt

This process will guide you in creating a new empty local archive database, if the current database has been corrupted. A typical symptom of a corrupt Local Archive database is that you are not able to connect to **LocalArchive-Int.HD** dataflow.

- Use service dongle to log onto Windows on scanner.
- Stop the SQL Anywhere – LocalArchive10 service
- Move or rename the database files (.db, .dbs, .txt) from E:\Archive\LocalArchive (Images and reports are stored within Files folder)
- Restart system. An empty database will now be created.

**Note!** Ensure that the Patient Archive is backed up regularly to avoid loss of data