

PC Setup Instructions

This document details the pc work station requirements typically used with Clinician. This setup should be completed on all computers that create orders or print results.

System Requirements:

- High speed Internet
- Internet Explorer 10 or above is preferred (older browsers are backward compatible however, they are no longer supported by Microsoft)
- Adobe Acrobat Reader is required for PDF results.

The next 2 steps **require ADMIN rights** on the PC for installation. Once installed, any user will have access to them.

1. Download the latest version of Meadco Script X directly from the Meadco website [HTTP://WWW.MEADROID.COM/SCRIPTX/USERLINK.ASP](http://www.meadroid.com/scriptx/userlink.asp)
OR click on the Download ScriptX here hyperlink that can be found on the top right of the Clinician Home page <https://clinician.changehealthcare.com>
2. Once installed, verify the version by clicking this link:
<http://scriptxsamples.meadroid.com/Basic/InstalledVersion/release/>
3. When prompted, accept the DX Control. The Dx Control will typically show up the first time a requisition is being printed. There is no external link to download this control.



Printer Settings (printer settings are found in the browser)

The steps below explain how to change the page setup in IE.


Log into **IE**, click on **File** in the upper left corner of the Internet Explorer Browser.

Click on **Page Setup**.

Add a check mark to Print Background Colors and Images

Set all four margins to .25".



Changing the Headers and Footers to empty is recommended.






Additional settings:

1. Under **TOOLS-Internet Options-Browsing History-Settings**, AUTOMATICALLY should be selected.






2. Click on **ADVANCED**- the items in yellow below should be Checked ✓. Click apply to save the settings.





Make sure these are also checked under the security section.





3. Compatibility View-typically only required when printing **RESULTS** directly from Clinician.

