

The screenshot shows the 'New Patient Registration' screen on an iPad. A left-hand sidebar contains a list of menu items: 'Demographics' (highlighted in teal), 'Responsible Party', 'Emergency Contact', 'Appointment Reminders', 'Photo', 'Insurance Cards', 'Allergies', 'Medications', 'General History', and 'Submit Form'. The main content area is titled 'Demographics' and includes a font size selector (A A A). Below the title is a paragraph: 'Please verify and update the below demographic information for the patient. If nothing has changed, please click the Confirm Section button.' The form contains several input fields: 'First Name*' (Kayla), 'Middle Name' (Brooke), 'Last Name*' (Jones), and 'Suffix'. There is also a 'Nick Name' field. The 'Street Address*' field contains '8169 Ashland Drive', and the 'Apartment/Suite #' field is empty. A 'Scroll for More' button with a downward arrow is located below the apartment field. At the bottom of the screen are two buttons: '< Previous Section' and 'Confirm Section >'. The iPad is shown in a black frame with a home button on the right side.

New Patient Registration

Font Size: A A A

Demographics

Please verify and update the below demographic information for the patient. If nothing has changed, please click the Confirm Section button.

First Name* Middle Name Last Name* Suffix

Kayla Brooke Jones

Nick Name

Street Address*

8169 Ashland Drive

Apartment/Suite #

Scroll for More ↓


< Previous Section Confirm Section >

MyMedicalLocker Intelligent Intake

iPad Setup Guide

Downloading the Application

1. Open up the App Store on the iPad
2. Tap the search button in the lower left hand corner.
3. In the search box type in “Intelligent Intake”
4. Select the one labeled “Intelligent Intake”. From there select “Get” or “Install”




Setup Application Settings

1. Navigate to the iPad settings
2. Scroll down to the section labeled “MyMedicalLocker Intelligent Intake”

ID: This setting allows you to give the iPad a specific identifier that can be used on the “Intake Scheduler” window.


Params: This setting is only for development purposes. This value needs to remain **empty**.




Setting up Guided Access

Guided Access allows you to limit your iPad to a single app and lets you control what features are available. We highly encourage using this method to “lockdown” iPads when using the MyMedicalLocker Intelligent Intake.


1. From the iPad settings navigate to the section labeled “**General**”. You should see new options on the right appear. Select the section labeled “**Accessibility**”.



2. Once inside “**Accessibility**” you should see quite a few different settings. The one we are going to activate is called “**Guided Access**”; turn this to “**On**”.



3. Tap the passcode settings, the tap “**Set Guided Access Passcode**”. Enter a passcode and then reenter it.



4. To start a guided access session. Open up the MyMedicalLocker Intelligent Intake app, then triple click the Home button. Then tap start to begin the session.

