

# Zedmed Cloud: Pathology and imaging providers

**This document provides information for pathology and radiology providers that need to setup laboratory software (e.g., Sonic) for a practice on Zedmed Cloud.**

Zedmed Cloud customers use a supported web browser to access Zedmed's practice management software. To facilitate secure messaging and files exchanges between Zedmed Cloud and external systems, middleware is installed on a local practice workstation.

When setting up secure file transfer software, the results files (HL7, ORU, etc) should be saved on the local machine running Zedmed's middleware.

The location the files need to be saved to is ***C:\ZedmedMiddleware\Results***

Files in this folder will then be uploaded to Zedmed by the middleware.